



Centrifuge Low Speed

- Double electronic lock.
- Memory program.
- Temperature protection & overshock protection.

Technical Data

Model	CLS-ELX-400
Rotor Capacity	6X15/50ml Centrifuge Tube 6X2/5/10ml Blood Tube
Speed	500 - 4000rpm
Max RCF	2146xg
Semidiameter Of Rotor	12cm
Timing Range	30sec - 59min 50sec
Power	80W
Power Supply	AC 110V - 220V, 50/60Hz
Dimensions (mm)	170 X 260 X 110 (WxDxH)