

Low Temp. Dry Bath Shaking Incubator



- Hot cover heat preservation function.
- Multi-point operation.
- Heating & cooling.
- Multi modules optional.

Technical Data

Model	LTDB-ELX-15
Actual temp range	RT AMBIENT TO -15°C ~100°C
Temperature Stability	±0.5°C
Max heating rate	5.5°C /min
Max cooling rate	5°C/min (100°C-Room Temp) 1°C/min (Below Room Temp)
Program memory	P1 - 9
Gradient Run	Yes
Speed	100-1500 rpm (Low Module)
Timing	5sec-24hrs
Power	150W
Voltage	100-240V 50/60HZ
Dimensions (mm)	265 X 190 X 150 (WxDxH)

Dry Bath Adaptors

Adaptor Options	Adaptor Capacity
ELX - 0.2	96X0.2ml Centrifuge Tube / 96-hole PCR Microplate / PCR Strips 8x0.2ml)
ELX - 0.5	24X0.5ml Micro Centrifuge Tube
ELX - 1.5	24X1.5ml Micro Centrifuge Tube
ELX - 2.0	24X2.0ml Micro Centrifuge Tube
ELX - 15	8X15ml Centrifuge Tube
ELX - 50	4X50ml Centrifuge Tube
ELX - Thermal	Thermal Lid Optional For ELX - 0.5 / 1.5 / 2.0ml