



Mini Centrifuge High Speed

- Enclosed high-speed type rotor.
- Stable rotation.
- Temperature protection function.
- Memory program.

Technical Data

| | |
|------------------------------------|--|
| Model | MCHS-ELX-145 |
| Control Mode | Digital Display & Button |
| Max Relative Centrifugal Force | 14104xg |
| Speed Range | 1000-14500rpm \pm 5% |
| Rotor Capacity | 0.2/0.5/1.5/2ml*12 Centrifuge Tube 0.5/1.0ml*12 Micro Blood Tube |
| Max Speed Acceleration / Reduction | \leq 12sec / \leq 15sec |
| Noise Under Max Rotation Speed | \leq 47db |
| Input Power | 55W |
| Shell Protection Level | AC100-240V / DC24V |
| Allow Surrounding Temperature | 5 - 40°C/ 80%RH |
| Dimensions (mm) | 280 X 390 X 220 (WxDxH) |