

Mini Dry Bath



- Small volume.
- Multi-point operation.
- LCD display.
- Over temperature protection.

Technical Data

| | |
|------------------------------|-------------------------|
| Model | MDB-ELX-16 |
| Temperature Range | Room temp + 5°C ~ 100°C |
| Temperature Stability | ±0.5°C |
| Temperature Homogeneity | ±0.5°C |
| Display Accuracy | 0.1°C |
| Heating Time | ≤12min (25-100°C) |
| Temperature Point Can Be Set | 9 |
| Time Setting Range | 5sec - 23hrs 59min |
| Power | 48W |
| Voltage | DC 24V |
| Dimensions (mm) | 160 X 115 X 96 (WxDxH) |

Mini Dry Bath Adaptors

| | |
|-----------------|-----------------------------------|
| Adaptor Options | Numbers And Aperture Of Test Tube |
| ELX-0.2 | 40 X 0.2ml PCR Tube |
| ELX-0.5 | 24 X 0.5ml PCR Tube |
| ELX-1.5 | 15 X 1.5ml Micro Centrifuge Tube |
| ELX-2.0 | 15 X 2.0ml Micro Centrifuge Tube |