



## Multi Tube Vortex Mixer

- Timing function.
- Large handling capacity.

### Technical Data

Model	MTVM-ELX-10
Motor Type	Brushless Motor
Horizontal Gyration Amplitude	4mm
Speed	1500 - 3000rpm
RPM Increment	10rpm
Timing	1min - 8hrs
Adaptors	(ELX101) Centrifuge tube : 25X50ml   16X15ml (ELX102) Centrifuge tube : 49X15ml   49X10ml (ELX103) Blood tube : 49X2/5ml (ELX104) Centrifuge tube : 81X1.5/2ml
Allowed Ambient Humidity Temp.	5 ~ 40°C 80%RH
Voltage	AC 100 - 220V, 50/60Hz
Power	30W
Dimensions (mm)	355 X 260 X 145 (WxDxH)