

## Multifunction Vortex Mixer



- Touch / Run, easy to switch.
- Various fixtures, large handling capacity.
- Brushless motor with strong power.
- Ultra-low body, stable operation.

### Technical Data

Model	VM-ELX-11
Operating Mode	Touch / Run
Control Mode	Digital Display & Button
Display	LCD
Motor	Brushless Motor
Oscillation	4.5mm (Circular Oscillation)
Speed Range	Touch: 0-3000rpm   Run: 100-3000rpm
Timing / Time	1min - 8hrs
Shell Protection Level	IP21
Voltage	20W
Power	DC12V
Dimensions (mm)	150 X 150 X 100 (WxDxH)

### Multifunctional Mixer Adapted Module

ELX01			Centrifuge tube / Test tube
ELX02			Flask / Multi tube
ELX03			21X12mm Blood collection tube
ELX04			12X5ml/15ml Centrifuge tube
ELX05			4X50ml Centrifuge tube   5X5ml/ 15ml Centrifuge tube
ELX06			96X0.2ml Centrifuge tube   96 well plate
ELX07			20X0.5ml Centrifuge tube   30X1.5ml/2.0ml Centrifuge tube
ELX08			Cell culture plate