



## Orbital Shaker

- Brushless motor.
- Suitable for a variety of containers.

### Technical Data

Model	OS-ELX-12	Timing	1min - 23hrs 59min
Type	Orbital	Max Load	2kg
Plate Size	250X250 mm	Allowed Ambient Humidity Temp.	5 ~ 40°C 80%RH
Rotation Angle	Circular	Voltage	AC 100 - 240V
Orbit Diameter	20	Power	30W
Speed	40-200 RPM	Dimensions (mm)	355 X 260 X 100 (WxDxH)