



3D Shaker

- Timing function.
- Suitable for a variety of containers.
- Adjustable speed.
- 3D mixing.

Technical Data

Model	3DS-ELX-13	Motor Type	Brushless Motor
Type	3D Oscillating	Timing	1min - 8hrs
Display	LCD	Allowed Ambient Humidity Temp.	5 ~ 40°C 80%RH
Tilt Angle	10°	Voltage	DC 12V
Max Load	0.8kg	Power	10W
Speed	10 - 80RPM	Dimensions (mm)	184 X 184 X 130 (WxDxH)